

## PAPER INSERTS

<u>Item #</u>	<u>Item Name</u>	<u>Aperture</u>	<u>List Price</u>
<a href="#"><u>0012-2479</u></a>	Paper Inserts With 1.5 mm Dia. Aperture - pk/200	Dia. 1.5 mm	\$215.60
<a href="#"><u>0012-2480</u></a>	Paper Inserts With 3.0 mm Dia. Aperture -pk/200	Dia. 3.0 mm	\$215.60
<a href="#"><u>0012-2481</u></a>	Paper Inserts With 5.0 mm Dia. Aperture - pk/200	Dia. 5.0 mm	\$145.00
<a href="#"><u>0012-2482</u></a>	Paper Inserts With 11.0 mm Dia. Aperture -pk/200	Dia. 11.0 mm	\$145.00
<a href="#"><u>0012-2483</u></a>	Paper Inserts With 1.5mm x 9mm Rectangular Aperture -pk/200	1.5 x 9 mm	\$160.00

Use of paper inserts facilitates mounting small samples and pressing KBr pellets. Paper inserts are used to make KBr pellets by dropping them into the bore of a 13mm die set. These inserts frame the pellet for ease of handling and also provide a visual reference for the proper amount of KBr powder which improves pellet quality and consistency.

### **INTERNATIONAL CRYSTAL LABORATORIES**

11 Erie Street, Garfield, NJ

[www.internationalcrystal.net](http://www.internationalcrystal.net)

[icmail@internationalcrystal.net](mailto:icmail@internationalcrystal.net)

Phone: 973-478-8944.

Fax: 973-478-4201

