

## Pres-N-Lok Demountable Cell

<u>Item #</u>	<u>Item Name</u>	<u>List Price</u>
<u>0006-2467</u>	Pres-N-Lok Cell Complete with Assorted Spacers, but without Windows	\$181.00
<u>0006-4407</u>	Pres-N-Lok Cell O-Rings-25mm, pk/12	\$59.50

Designed for optimum performance, ICL's Pres-N-Lok cell enables the user to analyze mulls and smears without the use of clumsy springs, nuts, bolts etc. Assembly of the cell requires only two crystal windows 25mm diameter by 4mm thick, one PTFE spacer and the user's sample.

### INTERNATIONAL CRYSTAL LABORATORIES

11 Erie Street, Garfield, NJ. 07026

[www.internationalcrystal.net](http://www.internationalcrystal.net)

[iclmail@internationalcrystal.net](mailto:iclmail@internationalcrystal.net)

Phone: 973-478-8944.

