

Real Crystal® IR Sample Cards

<u>Item #</u>	<u>Item Name</u>	<u>List Price</u>
<u>0000-7095</u>	Real Crystal® IR Cards- NaCl Pack of 25, 9.5mm aperture	\$138.00
<u>0000-7096</u>	Real Crystal® IR Cards- KBr Pack of 25, 9.5mm aperture	\$170.00
<u>0000-7097</u>	Real Crystal® IR Cards- KCl Pack of 25, 9.5mm aperture	\$147.00
<u>0000-7131</u>	Real Crystal® IR Cards- NaCl Pack of 25, 15mm aperture	\$182.00
<u>0000-7132</u>	Real Crystal® IR Cards- KBr Pack of 25, 15mm aperture	\$182.00
<u>0000-7133</u>	Real Crystal® IR Cards- KCl Pack of 25, 15mm aperture	\$185.00
<u>0000-7087</u>	Real Crystal® IR Cards- NaCl Pack of 25, 19mm aperture	\$170.00
<u>0000-7090</u>	Real Crystal® IR Cards- KBr Pack of 25, 19mm aperture	\$201.00
<u>0000-7091</u>	Real Crystal® IR Cards- KCl Pack of 25, 19mm aperture	\$177.00
<u>0000-7092</u>	Real Crystal® IR Cards- NaCl Cover Slips, Pack of 25	\$118.00
<u>0000-7093</u>	Real Crystal® IR Cards- KBr Cover Slips, Pack of 25	\$138.00
<u>0000-7094</u>	Real Crystal® IR Cards- KCl Cover Slips, Pack of 25	\$126.00
<u>0000-7086</u>	Real Crystal® IR Cards- NaCl Sample Pack of 5, 9.5mm aperture	\$73.00
<u>0000-7088</u>	Real Crystal® IR Cards- KBr Sample Pack of 5, 9.5mm aperture	\$73.00
<u>0000-7089</u>	Real Crystal® IR Cards- KCl Sample Pack of 5, 9.5mm aperture	\$73.00

<u>Item #</u>	<u>Item Name</u>	<u>List Price</u>
<u>0000-7134</u>	Real Crystal® IR Cards- NaCl Sample Pack of 5, 15mm aperture	\$73.00
<u>0000-7135</u>	Real Crystal® IR Cards- KBr Sample Pack of 5, 15mm aperture	\$73.00
<u>0000-7136</u>	Real Crystal® IR Cards- KCl Sample Pack of 5, 15mm aperture	\$73.00
<u>0000-7098</u>	Real Crystal® IR Cards- NaCl Sample Pack of 5, 19mm aperture	\$73.00
<u>0000-7099</u>	Real Crystal® IR Cards- KBr Sample Pack of 5, 19mm aperture	\$73.00
<u>0000-7100</u>	Real Crystal® IR Cards- KCl Sample Pack of 5, 19mm aperture	\$73.00

ICL's new Real Crystal® IR sample cards provide a conveniently mounted non-absorbing substrate for qualitative infrared analysis. They are ideal for liquids, mulls, film casting from solutions and pyrolysis samples. The 2" x 4" cards are available with 9.5mm, 15mm and 19mm clear apertures and they will fit in the universal slide mount of any spectrophotometer. Space is provided to write on the card to identify the date of the analysis, the sample identity and the name of the analyst.

INTERNATIONAL CRYSTAL LABORATORIES

11 Erie Street, Garfield, NJ
www.internationalcrystal.net
iclmail@internationalcrystal.net
 Phone: 973-478-8944.
 Fax: 973-478-4201

